



# Brachytherapy sources



<p>DOSEO platform includes two types of brachytherapy sources: <math>^{125}\text{I}</math> seeds and a <math>^{192}\text{Ir}</math> projector.</p>		
	$^{125}\text{I}$ seeds	$^{192}\text{Ir}$ projector
Average energy	20 keV	380 keV
Initial maximum activity	740 MBq	370 GBq
Reference air kerma rate	$30 \mu\text{Gy}\cdot\text{h}^{-1}$ at 1 centimeter	$30 \text{mGy}\cdot\text{h}^{-1}$ at 1 meter